

# Instructions for Users

## Model: PULSAR PRO HIGH

### PULSAR PRO HIGH

PULSAR PRO HIGH is intended for the support of mild to moderate ankle distortions, Tibialis Tendinitis, Osteoarthritis of ankle and foot, Chronic ankle instability, Tendinopathy of foot and ankle; Prophylactic and therapeutic use in sports activities. Rupture- tendon of the foot and ankle, non-traumatic. Acquired flat-foot or pronation of the foot/ankle.



1. Open all the straps and slide the injured ankle into the brace.



2. Put the foot of the injured ankle upon the bottom foot plate and secure the strap on the foot.



3. Locate the middle strap and feed through the D-Ring on the opposite side and attach to desired support. The length of this strap is adjustable by shortening or lengthening the strap on the underside.



4. Locate the top strap and feed through the D-Ring on the opposite side and attach to desired support. The length of this strap is also adjustable in the same way as the middle one. DO NOT OVERTIGHTEN.



### CARE AND CLEANING INSTRUCTIONS

Hand wash with a mild soap and let air dry.

### PRO ANKLE BRACE SIZE CHART

| Size       | Women's Shoe Size | Men's Shoe Size |
|------------|-------------------|-----------------|
| S/M Left   | Up to Size 10     | Up to Size 9    |
| S/M Right  | Up to Size 10     | Up to Size 9    |
| L/XL Left  | 10.5+             | 9.5+            |
| L/XL Right | 10.5+             | 9.5+            |