

SIZE INFORMATION

Thigh Sizing Guide- measure 15cm (6") above the medial condylar plateau.
Calf Sizing Guide- measure 12cm (5") below the medial condylar plateau.

Size	Thigh Circum. (CM)	Calf Circum. (CM)
S	33-39.5	27-31
M	39.5-47	31-35
L	47-53.5	33.5-38
XL	53.5-59.5	36-40

WASHING INSTRUCTION

- Remove metal stays and close all straps before washing.
 - Hand wash in cold water with mild soap.
 - Dry clean, bleach or tumble is prohibited.
 - Natural dry, and do not expose this brace to any heat source.
- Occasionally apply a dry lubricant on hinge especially after expose to water.



STORAGE CONDITION

The product should be stored in a room with a temperature of -5°C to 40°C, a humidity of 35% to 75%, good ventilation, dry, clean, and avoid direct sunlight.

INDICATION

- Knee instability due to anterior and/or posterior cruciate ligament (ACL/PCL) deficient
- To restrict the range of motion after reconstructive surgery
- Recommended for low to medium impact sports

INTENDED USE

To offer immobilization and/or controlled range of motion rehabilitation to the knee joint post surgery.

CONTRAINDICTION

- Paraesthesia
- Severe varicose vein
- Skin disorders or injuries in the relevant area
- Impaired lymphatic drainage

CAUTION

Please consult with your physician if any allergic symptom occurs. Do not use more than 24 hours of continuous.

RECOMMENDATION

These products are to be used only by those which have been medical specialist. In order to dispose of the product, comply with the legal norms of your country. If any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

CLASSIFICATION

Class I medical device, non sterile, no measurement function.
This product fulfils the Medical Device Regulation (EU) 2017/745.




Application Instruction

KNEE BRACE



INSTRUCTION MANUAL


KNEE BRACE




1. Measure thigh and calf circumference with leg extended. From sizing chart to select proper size.

2. Unfasten all the straps. Sit on the ground with your knees bent at 30 degrees. Adjust all the inner linings to be centered on the support frame. Put on the knee pads and pull them down to the knees.


3. Place hinge axis in line up with knee rotational axis.



4. Press the rigid shin shell against the leg and fasten strap 3 behind the hamstring. Next, fasten strap 2 behind the calf.





5. Be sure that the anterior thigh shell is pressed flat against the leg and fasten strap 4 in the back above knee joint. Next, fasten anterior calf strap 5 in the front below the knee; then fasten posterior calf strap 1 in the back below knee. Be sure to fasten all of the straps in the proper order.



Note: It is important to fasten the straps in the order described above. After 10-15 minutes of activity, re-adjust all the straps in the same order as described. This is required because of tissue swelling and compression.

How to Set the Hinge



1. Use a screwdriver to remove screws for inserting flexion/ extension stops.
2. Take off the existing stop and insert new plug-in stops, which are desired flexion/extension degree.
3. Screw it back to ensure settings. Be sure both screws are seated and snug.